

Waste Amalgam Storage

The PHS Wastemanagement Waste Amalgam container is designed to help protect staff from the potentially harmful effects of mercury inhalation, by containing an integral mercury suppressant foam.

Product Information	
Capacity	650ml
Size	105mm x 100mm (dia)
EWC Code	18 01 10

Instructions for use:

Simply store amalgam waste in this container until ready for collection by your PHS representative.



Amalgam Capsule Storage

PHS Wastemanagement offers a range of robust, sealed containers to allow the secure storage of dental amalgam capsule waste, with an integral mercury suppressant to combat the risk of vapour inhalation.

Product Information	
Capacity	1.8L or 6L
Size	170mm x 190mm (dia) 260mm x 200mm (dia)
EWC Code	18 01 10

Instructions for use:

Simply store used amalgam capsules in this container until ready for collection by your PHS representative.