

# Cytotoxic Waste Containers

Helping your premises stay safe and complaint

**phs** provides a range of specialist containers to allow you to effectively identify and segregate any cytotoxic or cytostatic waste you may produce. These units are highly visible and clearly labelled to encourage the correct identification and segregation.



## Instructions for use

These containers should only be used for cytotoxic and cytostatic hazardous waste, ensuring that items are safely contained within the unit. When ready for disposal, ensure the lid is securely sealed.

## Compliant and secure

The highly visible yellow units are constructed from durable and hard-wearing polypropylene, are puncture, impact, fluid and solvent resistant and approved to UN3291. These units are designed for easy use and closure once ready for disposal. Users should follow instructions and always treat cytotoxic waste carefully, keeping it segregated from other waste. EWC code 18 01 08.

## Key features:

- Approved to UN3291
- Puncture, impact, fluid and solvent resistant
- Highly visible yellow unit with purple lid in accordance with legislation

# Sizing

| Size (Ltr) | Capacity to fill line (Ltr) | Weight (g) | Dimensions (mm) |       |        | Aperture (mm) |       | Label |
|------------|-----------------------------|------------|-----------------|-------|--------|---------------|-------|-------|
|            |                             |            | Height          | Depth | Length | Length        | Width |       |
| 1          | 0.85                        | 126        | 187             | 108   | 90     | 58            | 36    | IML   |
| 3          | 2.21                        | 250        | 203             | 126   | 197    | 110           | 55    | IML   |
| 4          | 3.20                        | 320        | 245             | 175   | 175    | 110           | 60    | IML   |
| 13         | 10.50                       | 709        | 353             | 292   | 195    | 210           | 65    | IML   |
| 24         | 18.85                       | 1073       | 450             | 230   | 335    | 365           | 95    | IML   |

For further information please contact us on :

United Kingdom : 029 2080 9098 | [productinfo@phs.co.uk](mailto:productinfo@phs.co.uk) | [phs.co.uk](http://phs.co.uk)

Republic of Ireland : 01 643 4680 | [productinfo@phswashrooms.ie](mailto:productinfo@phswashrooms.ie) | [phswashrooms.ie](http://phswashrooms.ie)

